
Instructions

Self-testing kit for the diagnosis of urinary tract infections

SPECIFICATIONS

The Geratherm® URINARY TRACT INFECTION TEST consists of a test strip with 4 reaction fields. These can quickly check for the following parameters: Leucocytes, nitrite, blood and protein in the urine.

REAGENTS AND DETECTION LIMITS

Reagent	Test duration	Composition	Detection limits in urine
Leukocytes (LEU)	2 minutes	Amino acid ester pyrazole derivative, diazonium salt, buffer solution, non-reactive ingredients.	9 leukocytes/µl
Blood (BLO)	1 minute	3,3',5,5'-tetramethylbenzidine (TMB), diisopropylbenzole dihydroperoxide, buffer solution, non-reactive ingredients.	Haemoglobin: 0.018 mg/dl or 5 Ery/µl (with ascorbic acid < 50 mg/dl)
Nitrite (NIT)	1 minute	P-Arsanilic acid, N-1-(1-naphthyl)-ethylene-diamine, non-reactive ingredients.	0.05 mg/dl (with ascorbic acid < 30 mg/dl)
Protein (PRO)	1 minute	Tetrapromphenol blue, buffer solution, non-reactive ingredients.	For albumin: 0.075 g/l

PRECAUTIONS

Please read the instructions carefully before conducting the test.

- For home use.
- Only to be used as a urine test, not as a blood test.
- Do not use after the expiry date.
- Keep out of children's reach.
- For in-vitro diagnosis, not for internal application.
- After use, the test kit must be disposed of in compliance with local waste disposal regulations.

STORAGE AND USE

- Store in a dry place at a temperature between 2 °C and 30 °C.
- Do not freeze; do not expose to direct sunlight.
- Keep the test strips in the sealed pouch until ready to use.
- If the test strips are exposed to the air for a long period, the reagents may be destroyed, leading to false test results.
- Do not use the test strips if the pouch is torn or damaged.
- Please do not touch the reaction fields on the test strips.
- If the test strips are faded, discard them; they may be defective.

SCOPE OF DELIVERY


Materials provided	• Colour chart	• Instructions
• Urine test strips	Materials required but not provided	
• Sample collection container	• Clock or timer	• Kitchen paper

INSTRUCTIONS FOR USE


Before testing ensure that the test strip and the urine sample are at room temperature (15 °C - 30 °C). It is best to use the first urine of the morning and to take the sample midstream.

- Please take the necessary hygiene measures before conducting the test.
- Do not remove the test strip from the aluminium pouch until just before use, and grasp it as far away from the reaction fields as possible.
 - Do not hold the test strip in the urine stream until it has been flowing for 1-2 seconds. Ensure that all the reaction fields are thoroughly immersed in the urine (see Fig. 1).

Note: You can also collect the urine in a clean, dry container. In this case, dip the reaction fields completely in the urine and remove the test strip again immediately in order to prevent the reagents dissolving. When removing the strip from the container, carefully wipe it on the edge of the container to remove excess urine (see Fig. 1).

- Hold the test strip horizontal and lay it on a piece of absorbent tissue (e.g. a paper handkerchief) in order to prevent the chemical substances on the reaction fields and/or dirty hands coming into contact with the urine (see Fig. 2).
- Carefully compare the reaction fields with the corresponding colour bars on the colour chart. To do this, hold the test strip next to the colour bars (see Fig. 3).

ASSESSMENT		نتائج			
Reading of coloured test strips	Reagents	Results			
	Leucocytes (LEU)	- or ±	- or ±	+, ++ or +++	+, ++ or +++
Nitrite (NIT)	-	+	+	-	
Result	Leucocytes: Negative Nitrite: Negative	Leucocytes: Negative Nitrite: Positive	Leucocytes: Positive Nitrite: Positive	Leucocytes: Positive Nitrite: Negative	
Assessment	You probably do not have a urinary tract infection.	You probably have a urinary tract infection.	You probably have a urinary tract infection.	Carry out the procedure again with a new strip after a thorough personal cleaning.	
Recommendations	Should the symptoms continue, consult a doctor.	Consult a doctor.	Consult a doctor.	If you get the same results, consult a doctor.	

Certain medications or physiological states can lead to false positive or negative results. Never make a medical decision without consulting a doctor.

PRINCIPLES AND CRITICAL VALUES

LEUCOCYTES: The test procedure for leucocytes, or white blood cells, consists in the detection of the enzymes specific to these cells. These cells are present in an inflammation caused by an infection. The colour of the reagent zone changes from pink to violet if there is an infection. A raised protein concentration (> 500 mg/dl), a raised glucose concentration (> 2 g/dl) or a raised antibiotic concentration can lead to false results.

NITRITES: Nitrites are formed by the bacteria in question in most urinary tract infections (enterobacteriaceae). A urinary tract infection caused by these bacteria causes the reagent zone to become pink. Fasting, a vegetable-free diet or the intake of drugs containing phenazopyridine can lead to false results.

PROTEINS: Proteins in the urine can lead to false positive results in the LEUCOCYTES parameter. If there is protein in the urine, the reagent zone changes from yellow to green or also to blue.

BLOOD: If blood is present, the reagent zone in question changes from orange to green or also to dark blue. Blood in the urine is an indication of a serious problem. If the result is positive, consult a doctor. Tests carried out on women during their periods often show positive results.

PERFORMANCE DATA AND PROPERTIES

Accuracy: The samples were tested with a Geratherm® test strip for urinary tract infections. The same samples were tested for comparison with another product on the market (Bayer Multistix 10 SG Reagent Strips for Urinalysis).

Parameters	Results (Geratherm/Bayer)	Confidence interval of 95 %
Leucocytes	122/125	97,6 % (93 % - 99 %)
Nitrites	125/125	100 % (97 % - 100 %)
Protein	110/125	88 % (81 % - 93 %)
Blood	120/125	96 % (91 % - 99 %)

Practicability: A practicability test was carried out, in which the reading by patients was compared with the reading by trained professionals.

Reading of Geratherm® test strips for diagnosis of urinary tract infections by patients vs. trained staff.

Parameters	Results (Patient/Trained staff)	Confidence interval of 95 %
Leucocytes	100/100	100 % (96 % - 100 %)
Nitrites	100/100	100 % (96 % - 100 %)
Protein	100/100	100 % (96 % - 100 %)
Blood	100/100	100 % (96 % - 100 %)

BIBLIOGRAPHY

- Tietz NW. Clinical Guide to Laboratory Tests. W.B. Saunders Company, 1976.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry 2nd Ed. 2205, 1994.

LIST OF SYMBOLS

	Important: follow instructions		Number of test strips per test kit		Manufacturer
	In-vitro diagnostic agent		Expiry date		For single use only
	Min. and max. temperatures 2 °C - 30 °C		Lot number		reference number

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Geratherm Medical AG